



NCCU CLINICAL GUIDELINES  
SECTION: 4

THERMOGEGULATION

Section: 4 Thermoregulation  
References  
Date Revised: July 2011

Neonatology Clinical Guidelines  
King Edward Memorial/Princess Margaret Hospitals  
Perth Western Australia

---

## REFERENCES

---

- RPA ( Sydney)  
RWH ( Melbourne ) 2011  
Auckland Health 2010  
Great Ormond Street Hospital for Children review 2010  
ERBG guidelines UK 2006
- Agren, J., Sjors,G.,Sedin,G. (2006). Ambient Humidity Influences the Rate of Skin Barrier maturation in Extremely Preterm Infants. *Journal of Pediatrics* (148), 613-617.
- Bhatt, D. R., White,R., Martin,G., Van Marter,L.J.,  
Finer,F.,Goldsmith,J.P.,Ramos,R.,Kukreja,S.,Ramanathan,R. (2007). Transitional hypothermia in preterm newborns. *Journal of Perinatology* (27), S45-S47.
- Bohnhorst, B. (2010). Skin to Skin care in neonatal intensive Care Unit: more data regarding Seriously Ill Infants are badly Needed. *Neonatology* (97), 318-320.
- B;Lackburn et al (2001) Neonatal thermal care part III The effect of Infant position and temperature probe placement, *Neonatal network* 20/3/april 2001
- Cinar, N. D., Fils,T.M. (2006). Neonatal Thermoregulation. *Journal of Neonatal Nursing* (12), 69-74.
- Cramer, K., Wiebe, N., Hartling, L., Crumley, E., & Vohra, S. (2005). Heat Loss Prevention: A Systematic Review of Occlusive Skin Wrap for Premature Neonates. *Journal of Perinatology*, 25(12), 763-769.
- Darcy, A. E. (2009). Complications of the Late Preterm infant. *Journal of Perinatal Neonatal Nursing*, 23(1), 78-86.
- De La Fuente, L., Campbell,D.B.,Rios,A.,Grieg,A.,Graff,M.,Brion,L.P. (2006). Frequency analysis of ar and skin temperatures in neonates in servo controlled incubators *Journal of Perinatology*(26), 301 - 305.
- Elabbassi, E. B., Belghazi,K.,Delanaud,S.,LibertJP. (2004). Dry heat Loss in incubator: comparison of two premature newborn sized manikins. *Environment Toxique Perinatology - Adaptions Physiologiques*(92), 679-682.
- Healthcare, GE. Omnibed.
- Heimann, K., VaeBen,P.,Peschgens,T.,Stanzel,S.,Wenzel,T.G.,Orlikowsky.T.,. (2010). Impact of Skin to Skin Care, Prone Supine Positioning on Cardio respiratory Parameters and Thermoregulation in Premature Infants.

Hey, E. N., Katz,G. (1970). The Optimum Thermal Environment for Naked Babies. *Archives of Disease In Childhood*(45), 328.

Kim, S. M., Lee,E.Y.,Chen,J.,Ringer,S.A. (2010). Improved care and growth outcomes using hybrid humidified incubators in very preterm infants. *Pediatrics*, 125 (1): e 137-45.

Knobel, R. B., Holditch-Davis,D. (2006). Thermoregulation and heat loss prevention after birth and during neonatal intensive care unit stabilisation of extremely low -birthweight infants. *Advances in Neonatal Care*, 10(5), S7-S14.

Knobel, R. B., Holditch-Davis,D.,Schwartz.T.A. (2010). Optimal body temperature in transitional extremely low birth weight infants using heart rate and temperature as indicators. *JOGNN*, 39(3-14).

Knobel, R. B., Vohra, S., & Lehmann, C. U. (2005). Heat Loss Prevention in the Delivery Room for Preterm Infants: A National Survey of Newborn Intensive Care Units. *Journal of Perinatology*, 25(8), 514-518.

Laptook, A. R., Salhab,W., Bhaskar,B. . (2007). Admission temperature of low birth weight infants: predictors and associated morbidities. *Pediatrics*, 119(3), e643 -e649.

Laroia, N., Phelps,D.L.,Roy,J. (2007). Double wall versus single wall incubator for reducing heat loss in very low birth weight infants in incubators. *The Cochrane Collaborations* (4).

Lee, H. C., Ho, Q. T., & Rhine, W. D. (2008). A quality improvement project to improve admission temperatures in very low birth weight infants. *Journal of Perinatology*, 28(11), 754-758.

Libert, J., Bach,V.,Farges,G. (1997). Neutral temperature range in incubators: Performance of Equipment in current use and new developments. *Critical reviews in biomedical engineering*, 25, 287-370.

Mance, M., Short,M. (2008). Keeping infants warm: Challenges of hypothermia. *Advances in Neonatal Care* 8(1), 6-12.

McCall, E. M., Alderdice,F.,Halliday,H.L.,Jenkins,J.G.,Vohra,S. (2010). Interventions to prevent hypothermia at birth in preterm and/or low birthweight infant ( review). *Cochrane Data Base of systematic reviews* (3), 1-53.

Meyer, M. P., Payton,M.P.,Salmon,A.,Hutchinson,C.de Klerk,A. (2001). A clinical comparison of radiant warmer and Incubator care for preterm infants from birth to 1800gms. *Pediatrics*, 2(108), 395-401.

Meyer, M. P., Bold,G.T., . (2007). Admission temperatures following radiant warmer or incubator transport for preterm infants < 28 weeks: a randomised study. *Archives of Disease In Childhood: fetal and Neonatal* 92(4), F293 - F297.

Mok, Q. B., C.A.,Ducker,D.A.,McIntosh,N. (1991). Temperature instability during nursing procedures in preterm infants. *Archives of Disease In Childhood* (6), 783-786.

New, K., Flenady,V.,Davies,M.W. (2008). Transfer of preterm infants from incubator to open cot at lower versus higher body weight ratio. *The Cochrane Collaboration*(4).

New, K., Bogossian,F.East, C.,Davies,M.W. (2009). Practice variation in the transfer of premature infants from incubators to open cots in Australian and New Zealand neonatal nurseries: Results of an electronic survey. *International Journal of Nursing Studies*.

Sherman, T. I., Greenspan, J. S., St Clair, N., Touch, S. M., & Shaffer, T. H. (2006). Optimizing the neonatal thermal environment. *Neonatal Network: NN*, 25(4), 251-260.

Sinclair, L., Crisp, J., Sinn, J. (2009). Variability in incubator humidity practices in management of preterm infants. *Journal of Paediatrics and Child Health* (45), 535-540.

Sjörs, G., Hammarlund, K., & Sedin, G. (1997). Thermal balance in term and preterm newborn infants nursed in an incubator equipped with a radiant heat source. *Acta Paediatrica*, 86(4), 403-409.

Soll, R. F. (2008). Heat loss prevention in neonates. *Journal of Perinatology*, 28(S1), S57-S59.

Telliez, F., Bach, V., Delanaud, S., Leke, A., Abdiche, M., Chardon, K. (2001). Influence of incubator humidity on sleep and behaviour of neonates kept at stable body temperature. *Acta Paediatrica*(90), 998-1003.